

Instructions to install phenomics workflow on a new computer

Installation of *JupyterLab*

Note: If you are viewing these instructions in *JupyterLab*, this is already complete.

1. Install *Anaconda*
2. Install *JupyterLab* by the command:
`conda install -c conda-forge jupyterlab<\code>`
3. Install the bash kernel for *JupyterLab* by the command:
`conda install -c conda-forge bash_kernel`
4. Open *JupyterLab* using the command: `jupyter lab`

Installation of **GM0detector**

1. Clone `gmodetector_py` from Github using the command:
`git clone https://github.com/naglemi/gmodetector_py.git`

1. Create a conda environment that already has most dependencies using this command:
`conda create -q -n test-environment python=$TRAVIS_PYTHON_VERSION pytest codecov numpy pandas matplotlib tzlocal h5py pytables`
2. Activate the new conda environment
`source activate test-environment`
3. Install additional packages through conda:
Install *rpy* to allow loading of R packages in Python
`conda install -c conda-forge rpy2`
Install *GNU parallel* to allow parallelization
`conda install -c conda-forge parallel`
4. Install *Spectral Python* from the *pypi* repository
`pip install spectral`
5. Launch R and install *scales* from the R console via:
`install.packages('scales', repos = 'http://cran.us.r-project.org')`
6. Install *gmodetector_py* and test dependencies:
`pip install ".[test]"`

Installation of segmentation tools

Install all dependencies for deeplab in a new conda environment

```
conda create --name deeplab
```

```
conda activate deeplab
```

```
conda install -c anaconda scipy
```

```
conda install -c conda-forge opencv
```

```
conda install -c conda-forge tensorflow=1.14
```

```
conda install -c anaconda yaml # Not needed
```

```
pip install pyyaml
```

```
conda install -c anaconda pillow
```

Copy or clone the appropriate deeplab repo with pretrained model (not yet all in one place on Github at time of writing; currently in deeplab folder on OSU server steed)

Installation of tools for classifying missing/contaminated explants

Follow instructions on repo <https://github.com/Contamination-Classification/DenseNet> (<https://github.com/Contamination-Classification/DenseNet>) and note that tensorflow v1.14 can be used if there is an issue installing v1.13

Installation of image alignment

```
git clone https://github.com/NSF-Image-alignment/ImageAlignment.git
```

```
cd ImageAlignment
```

Create the environment and install dependencies available through pip

```
conda create -n alignment python=3.7
conda activate alignment
pip3 install pillow scikit-image numpy matplotlib easydict opencv-python pandas argparse xlrd
```

Install some additional dependencies through conda

```
conda install -y -c anaconda xlrd
conda install -y numpy
conda install -y -c conda-forge opencv
conda install -y -c conda-forge easydict
conda install -y scikit-image
conda install -y -c anaconda pandas
```

Installation of GM0labeler

1. Create and activate a new environment

```
conda create -n gmolabeler
conda activate gmolabeler
```

1. Install R. Instructions can be found here: <https://cloud.r-project.org/> (<https://cloud.r-project.org/>)
Alternatively, R may be installed through conda, either by installing R itself alone or by installing R as a dependency for rpy2 as shown below.

```
conda install -c conda-forge rpy2
```

1. Launch the R console and install required R packages from the R console:

```
install.packages("optparse")  
  
install.packages("data.table")  
  
install.packages("stringr")
```

1. Install Python (and R) dependencies

```
# Python dependencies
conda install -y -c anaconda pillow
conda install -y -c anaconda pandas
conda install -c anaconda pytables

# One R dependency installed by conda for randomcoloR R package
conda install -c conda-forge libv8
```

Install R package dependencies from within R

R

```
install.packages("scales")  
install.packages("readxl")  
install.packages("ggplot2")  
install.packages("randomcoloR")  
install.packages("tidyr")
```